



# Multifocal Contact Lenses with the Contex MF-19




**Inventors of the OK<sup>®</sup> Lens**

Copyright Contex, Inc. 1999-2007

**There are over 89  
Million Presbyopes  
in the U.S. alone!**



# Only 1% Currently Wear Contact Lenses



**Up to 40% of them  
May Require GP  
Lenses**



# These Numbers are Similar Worldwide




# Are You Ready to Fit these Patients?





# **The Simple-fit Multifocal-19<sup>®</sup>**

# Contex Multifocal-19<sup>®</sup>




**The “Comfortable” GP Multifocal**

Copyright Contex, Inc. 1999-2003



# Multifocal-19<sup>®</sup> Lens Design




# Multifocal-19<sup>®</sup> Advantages

- **Clear Distance Vision**
- **True “Aspheric” Add Power**
- **Highest Resolution Optics**
- **Comfortable Design**
- **Simultaneous Distance and Near Vision**


# Multifocal-19<sup>®</sup> Lens Design

## Simultaneous Vision



# **MF-19<sup>®</sup> Fitting Characteristics**

- **Slight Apical Clearance  
Progressing to Mid-Peripheral  
Touch**
- **Lens Must Position Central to  
Superior for Best Results**



**Select Patients  
with +2.00 Add  
Power or Lower**

# **Interview Patient Regarding:**

- **Lens Expectations**
- **Career and Hobbies**
- **Previous Lens History**

# Corneal Data Required

- **Central K Readings**
- **Refractive Error**
- **Eyelid Structure**

# 2 Fitting Methods

- **Central K's, Rx**
- **Trial Fitting**



# Central K's and Rx

- **Simple Chart Makes  
Selecting Proper Base  
Curve Easy**

# Central K's and Rx

- If ordering empirically with CK's and Rx the Base Curve should be  $2 \frac{3}{4}$  diopters steeper than the flat Central K
- Example: Flat CK = 43.00 (7.85)  
Base Curve = 45.75 (7.38)

# Refraction vs. Lens Power

- Determine the Non-Dominant eye and Over-Plus the Lens Power by  $+0.50$  for Maximum Add

Example: Rx =  $-3.00$

Calculated Distance Power  $-5.50$

Calculated Non-Dom Power  $-5.00$



# **Trial Fitting Offers Superior Results**


# Trial Fitting

- Perform spherical over-refraction over the best fitting lens
- Adjust for vertex if above  $\pm 4.00$  diopters
- Over refraction is best with loose trial lenses or a trial frame(not a phoropter)

# Look For

- **Central to slightly superior lens position to provide good distance vision and reduce flare and glare in dim light or at night**
- **A lens that translates easily across the vertical corneal surface as the patient looks from distance to reading tasks**


# Lens in Translation



# Fluorescein Pattern on Straight Gaze








**During downward  
gaze the lens should  
easily translate up to  
provide maximum  
reading vision**

# Desired Fluorescein Pattern

- **1-3 mm Central Clearance  
Progressing to Mid-Peripheral  
Touch**
- **.8 to 1.2 mm Wide Mid-Peripheral  
Bearing**
- **Adequate Edgelifit**
- **Lens Must Be Centered!**
- **1-3 mm Movement On The Blink**

# Typical MF-19<sup>®</sup> Fluorescein Pattern



# Troubleshooting Poor Distance Vision

- Decentered lens, can cause patient to look through intermediate zone
- **Solution:** Steepen Base Curve



# Troubleshooting Poor Near Vision

- De-centered and over-minusing due to patient looking through intermediate zone
- **Solution:** Steepen Base Curve to center lens and verify Over-Refraction



# Case Study #1

## Original Data

**R- 44.00/45.37      .48e/TK- 40.12**

**L- 44.37/45.62      .52e/TK- 40.25**

**Spec-Rx = -2.00 S**

**Spec-Rx = -2.25 S**

# Initial Right Lens

**MF-19<sup>®</sup>**

**7.22 / 10.0 / -4.50**

**Lens fit acceptable but riding a bit high. Steepen Base Curve**

# Second Right Lens

**MF-19<sup>®</sup>**

**7.14 / 10.0 / -5.00**

**Lens centers perfect and Add  
adequate**



# Right Eye Final Data

**Over Refraction  $-.25$**

**J7 on the reading card**

**Patient wears lenses**

**12 to 16 hours a day**

# Initial Left Lens

**MF-19<sup>®</sup>**

**7.14 / 10.0 / -4.25**

**Lens fit acceptable but riding a bit high. Steepen Base Curve**

# Left Eye Final Data

**Over Refraction  $-.50$**

**J7 on the reading card**

**Patient wears lenses**

**12 to 16 hours a day**

# Thank You For Your Attention



Questions  
&  
Answers

Copyright Contex, Inc. 1999-2007